

The VNL-J10 Scopes

High-definition and high performance at its best



The PENTAX Medical VNL-J10 product range offers a complete solution for advanced flexible endoscopy needs.

Used with PENTAX Medical HD video processors equipped with i-scan technology the new VNL series will provide support to both diagnostic as well as therapeutic procedures. The VNL product range can be used in a wide variety of applications in ENT such as FEES examinations and functional video stroboscopy examinations.

Get connected to
the world of HD



Advanced assessment

Better orientation and support of diagnosis due to

- latest chip technology
- improved design of distal tip with improved light distribution
- i-scan technology

Increased patient's comfort

- Smallest outer diameter for better tolerance
- Improved tapering of distal end for optimized insertion

Enhanced and safe diagnosis

High performance diagnostic & therapeutic devices

- New mechanical design of instrument channel for smooth insertion of endoscopic devices
- Bifurcated umbilical cable to allow for stroboscopy

Optimized patient and user safety

- Innovative single-use suction control valve design
- STERRAD compatibility

Easy endoscope control & excellent ergonomomy

- Newly designed control body
- Flexible light guide cable

DEFINA Scopes

	Routine VNL11-J10	Small VNL8-J10	Therapeutic VNL15-J10
Distal end diameter [mm]	3.5	2.4	4.8
Insertion tube diameter [mm]	3.6	2.9	4.9
Instrument channel diameter [mm]	–	–	2.0
Working length [mm]	300	300	300
Total length [mm]	500	500	575
Angle of view [°]	80	80	80
Focal range [mm]	3-50	3-50	3-50
Angulation up/down [°]	130/130	130/130	130/130

EMEA Headquarter
Germany



PENTAX Europe GmbH
Julius-Vosseler-Straße 104 . 22527 Hamburg
Tel.: +49 40 / 5 61 92 - 0 . Fax: +49 40 / 5 60 42 13
E-mail: info.emea@pentaxmedical.com . www.pentaxmedical.com

TÜV Süd CE0123 • Medical device class: IIa • This product must be used only by Healthcare professionals. Before usage and for detailed product specifications, please refer to the instructions for use. In the interest of technical process, specifications may change without notice.

PENTAX
MEDICAL
Excellence in Focus